

SPECIAL PROMO

VALID FROM MAY 1ST TO DECEMBER 31ST 2020


SOLID.LINE CLIMATE CHAMBERS AT A SPECIAL PRICE!

The new Solid.Line KBF-S climate chamber is a custom-made solution for long-term testing and accelerated stability testing for pharmaceutical products in accordance with the ICH Q1A guideline. The corrosion-free and condensation-free interior also ensures exceptionally stable temperature and moisture levels, even when fully loaded.

FEATURES:

- Temperature range: 0 °C to 70 °C
- Humidity range: 20% to 80% RH
- APT.line™ preheating chamber technology
- Humidity regulation with capacitive humidity sensor and steam humidification
- LCD to display temperature and humidity along with additional information and alarms
- Internal data logger, measured values can be read out in open format via USB

Series	115	240	720	1020
KBF-S	•	•	•	•



Model KBF-S 720

